



Neutrons
for
polymers

WMZ

Coh/Incoh
separation

SRO in
elastomer

How-it-works

QENS Con-
tamination

Larmor
Precession
Instruments

NSE

SESANS,
SERGIS

Dynamic
Nuclear
Polarization

Soft matter with polarised neutrons

Wojciech Zając



The Henryk Niewodniczański Institute of Nuclear Physics
Polish Academy of Sciences, Kraków

ESS Partner & Industry Day
Kraków, 25.03.2014



Residual strain in the PU-SiO₂ composite. WANS spectra (polarised neutrons) for PU.

- Neutrons for polymers
- WMZ
- Coh/Incoh separation
- SRO in elastomer
- How-it-works
- QENS Contamination
- Larmor Precession Instruments
- NSE
- SESANS, SERGIS
- Dynamic Nuclear Polarization





Neutron spin as a measuring probe

Neutrons
for
polymers

WMZ

Coh/Incoh
separation

SRO in
elastomer

How-it-works

QENS Con-
tamination

Larmor
Precession
Instruments

NSE

SESANS,
SERGIS

Dynamic
Nuclear
Polarization

- Polarized neutrons offer a unique opportunity of separating, **at the machine level**, coherent from incoherent scattering. See how it works:

- If $\hat{B} (\vec{I} \cdot \vec{\sigma})$ denotes the interaction operator between the neutron spin and the nuclear spin, then the following holds for the scattering intensity in the spin-flip and non-spin-flip channels:

$$I^{\uparrow\uparrow} = \frac{1}{3} |B|^2 I(I+1), \quad \text{and} \quad I^{\uparrow\downarrow} = \frac{2}{3} |B|^2 I(I+1)$$

Thus the outgoing scattering intensity from unpolarized nuclear spins is always one third without, and two thirds with flip of the neutron spin from the polarised beam. The coherent scattering I_{coh} and incoherent scattering I_{inc} are linearly related to the measured spin-flip and non-spin-flip scattering:

$$I_{corr}^{\uparrow\downarrow} = \frac{2}{3} I_{spin\ incoh} \quad \text{and} \quad I_{corr}^{\uparrow\uparrow} = I_{coh} + \frac{1}{3} I_{spin\ incoh} + I_{isotope\ incoh}$$

$$I_{coh} = I_{corr}^{\uparrow\uparrow} - \frac{1}{2} I_{corr}^{\uparrow\downarrow}, \quad \text{and} \quad I_{spin\ incoh} = \frac{3}{2} I_{corr}^{\uparrow\downarrow}$$

- Spin-incoherent scattering always implies neutron spin-flip, while coherent contains a combination of both spin-flip and non-spin-flip events.



Neutron spin as a measuring probe

Neutrons
for
polymers

WMZ

Coh/Incoh
separation

SRO in
elastomer

How-it-works

QENS Con-
tamination

Larmor

Precession
Instruments

NSE

SESANS,
SERGIS

Dynamic
Nuclear
Polarization

- Polarized neutrons offer a unique opportunity of separating, **at the machine level**, coherent from incoherent scattering. See how it works:
- If $\hat{B} \left(\vec{I} \cdot \vec{\sigma} \right)$ denotes the interaction operator between the neutron spin and the nuclear spin, then the following holds for the scattering intensity in the spin-flip and non-spin-flip channels:

$$I^{\uparrow\uparrow} = \frac{1}{3} |B|^2 I(I+1), \quad \text{and} \quad I^{\uparrow\downarrow} = \frac{2}{3} |B|^2 I(I+1)$$

Thus the outgoing scattering intensity from unpolarized nuclear spins is always one third without, and two thirds with flip of the neutron spin from the polarised beam. The coherent scattering I_{coh} and incoherent scattering I_{inc} are linearly related to the measured spin-flip and non-spin-flip scattering:

$$I_{corr}^{\uparrow\downarrow} = \frac{2}{3} I_{spin\ inc} \quad \text{and} \quad I_{corr}^{\uparrow\uparrow} = I_{coh} + \frac{1}{3} I_{spin\ inc} + I_{isotope\ inc}$$

$$I_{coh} = I_{corr}^{\uparrow\uparrow} - \frac{1}{2} I_{corr}^{\uparrow\downarrow}, \quad \text{and} \quad I_{spin\ inc} = \frac{3}{2} I_{corr}^{\uparrow\downarrow}$$

- Spin-incoherent scattering always implies neutron spin-flip, while coherent contains a combination of both spin-flip and non-spin-flip events.



The instrument

Neutrons
for
polymers

WMZ

Coh/Incoh
separation

SRO in
elastomer

How-it-works

QENS Con-
tamination


Larmor
Precession
Instruments

NSE

SESANS,
SERGIS


Dynamic
Nuclear
Polarization

Such opportunity offers the D7 Diffuse Scattering instrument(ILL):





The INS/QENS spectrum "anatomy"



Neutrons
for
polymers

WMZ

Coh/Incoh
separation

SRO in
elastomer

How-it-works

QENS Con-
tamination

Larmor
Precession
Instruments


NSE

SESANS,
SERGIS

Dynamic
Nuclear
Polarization



$$\text{EISF}(Q) = \frac{I_{el}(Q)}{I_{el}(Q) + I_{gel}(Q)} - \text{an example}$$





QENS on PMMA – an particularly nasty case of diffuse coherent scattering

Neutrons
for
polymers

WMZ

Coh/Incoh
separation

SRO in
elastomer

How-it-works


QENS Con-
tamination

Larmor
Precession
Instruments

NSE

SESANS,
SERGIS

Dynamic
Nuclear
Polarization





QENS on PMMA – an particularly nasty case of diffuse coherent scattering

Neutrons
for
polymers

WMZ

Coh/Incoh
separation

SRO in
elastomer

How-it-works


QENS Con-
tamination

Larmor
Precession
Instruments

NSE


SESANS,
SERGIS

Dynamic
Nuclear
Polarization





In NSE*, Neutron Spins Precess Before and After Scattering & a Polarization Echo is Obtained if Scattering is Elastic



* F. Mezei, Z. Physik, 255 (1972) 145





NSE principles

Neutrons
for
polymers

WMZ

Coh/Incoh
separation

SRO in
elastomer

How-it-works


QENS Con-
tamination

Larmor
Precession
Instruments

NSE

SESANS,
SERGIS

Dynamic
Nuclear
Polarization





The principle of **SESANS**

Neutrons
for
polymers

WMZ

Coh/Incoh
separation

SRO in
elastomer

How-it-works

QENS Con-
tamination


Larmor

Precession
Instruments

NSE


SESANS,
SERGIS

Dynamic
Nuclear
Polarization





Spin Echo Small Angle Measurement (SESAME) Spin Echo Resolved Grazing Incidence Scattering (SERGIS)



- 1: Specular reflection,
- 2, 3: Off-specular reflection (above horizon)
$$\theta = \theta_i + \sqrt{\theta_i^2 + 2n\lambda/d} \quad n = \pm 1$$
- 4: Off-specular diffraction (below horizon)
$$\theta = \theta_i + \sqrt{\theta_i^2 + 2n\lambda/d - Nb\lambda^2/\pi} \quad n = -1$$



Length scales probed by SESANS

Neutrons
for
polymers

WMZ

Coh/Incoh
separation

SRO in
elastomer

How-it-works


QENS Con-
tamination

Larmor
Precession
Instruments

NSE

SESANS,
SERGIS

Dynamic
Nuclear
Polarization



(Wei-Ren Chen et al., "The Theoretical Foundation of Spin-Echo Small-Angle Neutron Scattering (SESANS) Applied in Colloidal System", UCANS-II Indiana University; July 08th 2011; Bloomington, IN)



Principles of Dynamic Nuclear Polarization

Neutrons
for
polymers

WMZ

Coh/Incoh
separation

SRO in
elastomer

How-it-works

QENS Con-
tamination

Larmor
Precession
Instruments

NSE

SESANS,
SERGIS

Dynamic
Nuclear
Polarization

The scattering amplitude a of thermal neutron, carrying spin s ; ($|s| = \frac{1}{2}$) interacting with a nucleus of spin \mathbf{I} (Abragam, Goldman; 1982):

$$a = b + 2B (\mathbf{I} \cdot \mathbf{s})$$

The constants b and B are determined by the two eigenvalues $b_{(+)}$ and $b_{(-)}$ of the operator $(\mathbf{I} \cdot \mathbf{s})$, depending on whether \mathbf{I} and \mathbf{s} couple in the channels $I + \frac{1}{2}$ or $I - \frac{1}{2}$. They may be interpreted as **spin-independent** and **$\frac{1}{2}$ of spin-dependent** scattering lengths:

$$b = \frac{(I + 1) b_{(+)} + I b_{(-)}}{2I + 1}$$

$$B = \frac{b_{(+)} - b_{(-)}}{2I + 1}$$

and for the hydrogen atom they are: $b = -3.74$ fm, $B = 29.1$ fm.



Principles of Dynamic Nuclear Polarization

Neutrons
for
polymers

WMZ

Coh/Incoh
separation

SRO in
elastomer

How-it-works

QENS Con-
tamination


Larmor
Precession
Instruments

NSE

SESANS,
SERGIS

Dynamic
Nuclear
Polarization

We assume a polarized target with polarisation, P_N , along an axis and the neutron beam polarized along the same axis with polarisation P_n . Furthermore, there are no correlations between the relative orientations of two nuclear spins and their relative positions in an unpolarized target. Then the total polarization-dependent coherent and incoherent neutron scattering cross-sections are:



$$\sigma_{coh} = 4\pi [b^2 + P_n P_N \cdot 2IbB + P_N^2 I^2 B^2]$$

$$\sigma_{inc} = 4\pi [I(I+1)B^2 - P_n P_N I B^2 - P_N^2 I^2 B^2]$$

(J.H. Zhao et al, Physics Procedia 42 (2013) p.39-45)



Spin Contrast Variation

Neutrons
for
polymers

WMZ

Coh/Incoh
separation

SRO in
elastomer

How-it-works



QENS Con-
tamination

Larmor
Precession
Instruments

NSE

SESANS,
SERGIS

Dynamic
Nuclear
Polarization



Small-angle neutron scattering (SANS) intensity decreased for positive spin polarization and increased for negative spin polarization. (Noda et al., J. Appl. Cryst.

44 3 (2011) p.503-513)